

# Introduction to Liquid Chromatography

Columns  
System Components  
Applications  
Troubleshooting

Susan M. Steinike, M.S.  
HPLC Marketing Department  
February, 2006

# A Brief History of Chromatography

- 1903: Russian botanist Mikhail Tswett separates plant pigments
- 1938: Russian scientists Izmailov and Shraiber use “drop chromatography”, later perfected as Thin Layer Chromatography (TLC) by Kirchner in the U.S.
- 1952: Martin and Synge receive Nobel Prize for “invention of partition chromatography” or plate theory to describe column efficiency
- 1966: HPLC was first named by Horvath at Yale University but HPLC didn’t “catch on” until the 1970s
- 1978: W.C. Stills introduced “flash chromatography”, where solvent is forced through a packed column with positive pressure

# Modern HPLC


- Late 1970s/early 1980s
  - Instrumentation developed for high pressure solvent delivery: pumps, autosamplers, diode array detectors
  - More uniform packing material produced for columns
- Last 20 years
  - Nothing really “new”, but by returning to the basic theory of chromatography, even better columns are on the market: smaller particle sizes which yield faster separations, but require hardware to withstand higher pressures.

# What is Chromatography?

- Separation of a mixture into individual components.
- The separation uses a Column (stationary phase) and Solvent (mobile phase).
- The components are separated from each other based on differences in affinity for the mobile or stationary phase.
- The goal of the separation is to have the best **RESOLUTION** possible between components.



# The Most Basic Explanation of Chromatography Ever



# How Do You Get Separation?

- Hardware: pumps, injector, detector
- Column: particle diameter, column size, packing materials, and the dreaded equations
- Our seminar will focus on the contribution of each factor to perform separations.

# Outline

- Column Considerations
  - Theory (including, well...you know)
  - Different Stationary Phases
- Hardware Components
  - Pumps, Injectors, Detectors, etc.
  - Examples of Application-Specific Configurations
- Applications
  - Pharmaceuticals and Proteomics
  - Food and Beverage, Environmental
  - Research and Method Development


# Outline

- System Troubleshooting
  - Leaks, Reproducibility, Column Care, and More
- Chromatography Software
  - Method and Sequence Setup
  - Calibration Curves and Reporting
- Chromatography Hardware
  - Modular LC-20 Prominence
  - Integrated LC-2010HT

# Modern HPLC vs. Traditional LC Methods


- Classical open-column LC.
- Thin-Layer Chromatography (TLC) and paper chromatography.
- In modern HPLC the columns and packings are, in general, highly refined, high in resolving capacity, and are reusable.

# HPLC and Pre-HPLC Techniques



# Column Types

- Normal Phase LC
  - Polar stationary phase: Silica
  - Nonpolar mobile phase: Hexane, Ethyl acetate
  - The LEAST polar compound comes out first






# Normal Phase HPLC Columns

- Cyano: Rugged, moderate polarity, general use
- -OH (Diol): More polar and retentive
- Amino: Highly polar, less stable
- Silica: Very rugged, low cost, adsorbent (Unbonded)

The cyano column with a low polarity mobile phase (hydrocarbon with a small amount of another solvent) will act as a normal phase column.

# Column Types

- Reversed-Phase LC
  - Nonpolar stationary phase: C8, C18
  - Polar mobile phase: Water, ACN, Methanol
  - The MOST polar compound comes out first



# Reversed Phase HPLC Columns

- C-18, C-8: Rugged, general purpose, highly retentive
- C-3, C-4: Less retentive, used mostly for peptides & proteins
- Phenyl: Greater selectivity than alkyl-bonded
- Cyano: Moderate retention, normal & rev. phase
- Amino: Weak retention, good for carbohydrates

The cyano column with a high polarity mobile phase (Water/MeOH) will act as a reversed phase column.

# Normal vs. Reversed Phase

| <u>Parameter</u>          | <u>Normal Phase</u> | <u>Reverse Phase</u> |
|---------------------------|---------------------|----------------------|
| Polarity of Packing       | Medium to High      | Low to Medium        |
| Polarity of Solvent       | Low to Medium       | Medium to High       |
| Elution Sequence          | Low Polarity First  | High Polarity First  |
| Increase Solvent Polarity | Faster Elution      | Slower Elution       |

# Column Types

## ■ Ion Exchange LC

- Stationary phase contains charged groups
- SAX (Strong Anion Exchange):  $\text{NH}_3^+$
- WAX (Weak Anion Exchange):  $\text{NR}_2\text{H}^+$  (DEAE)
- SCX (Strong Cation Exchange):  $\text{SO}_3^-$
- WCX (Weak Cation Exchange): Carboxymethyl (CM)
- More highly charged analytes have stronger retention
- More “bulky” stationary phases have weaker retention

# Column Types


- Size Exclusion LC (also called Gel Permeation)
  - Stationary phase is a polymer (polystyrene-divinyl benzene or acrylamide) with a defined pore size
  - Large compounds cannot fit into the pores and elute first
  - Used to determine molecular weight distribution of polymers

# Typical Column Sizes

- Particle size: 5  $\mu\text{m}$ , 3  $\mu\text{m}$ , and smaller
  - Monodispersed means particles are the same size
  - Very important for stable pressure and flow
  - Smaller particles produce higher system pressure
- Pore size: 100-120  $\text{\AA}$  is typical
- Surface area: 300-350  $\text{m}^2/\text{g}$
- Carbon load: 9-12% for C8, 16-20% for C18
  - Higher carbon load = better resolution but longer run times
  - Lower carbon load = shorter run times, but may change selectivity vs. higher carbon load



# Idealized HPLC Separation



# Void Volume


- The void volume is the amount of “dead” volume in the column that is not taken up by the particles of stationary phase.
- In general, there is approximately 0.1 mL of void volume for each cm of column length, for columns with a 4.6 mm i.d. and 5  $\mu\text{m}$  particles

$$V_m \approx 0.5d_c^2L$$

Where  $V_m$  is the column volume in mL,  
 $L$  is the column length in cm, and  
 $d_c$  is the inner diameter in cm

# Void Volume

- The void volume is exactly determined by injecting a compound that is completely unretained, then using the chromatogram to calculate void volume.
- Elution time x flow rate = void volume



# What is Chromatography?

- Separation of a mixture into individual components.
- The separation uses a Column (stationary phase) and Solvent (mobile phase).
- The components are separated from each other based on differences in affinity for the mobile or stationary phase.
- The goal of the separation is to have the best **RESOLUTION** possible between components.

# Factors Influencing Resolution

Capacity Factor,  $k'$   
Selectivity Factor,  
 $\alpha$   
Efficiency,  $N$

# The Resolution Equation

- Resolution is defined as the completeness of separation from one analyte to another
- In general, resolution may be expressed as:

$$\begin{aligned}R_s &= 2(V_{rb} - V_{ra}) / (W_a + W_b) \\ &= 2(tr_b - tr_a) / (W_a + W_b)\end{aligned}$$

- Where  $V_{ra/b}$  = retention volume of peak a/b
- Where  $tr_{a/b}$  = retention time of peak a/b
- Where  $W_{a/b}$  = width of peak a/b

# Resolution

- For closely eluting or adjacent peaks, the resolution equation may be expressed as:

$$R_s = 1/4[(\alpha - 1) / \alpha] \sqrt{N} [k' / (1 + k')] ]$$

- The terms of capacity factor ( $k'$ ), selectivity ( $\alpha$ ), and efficiency ( $N$ ) all contribute to resolution
- Let's look at how each term affects resolution



# Capacity Factor, $k'$


The relative degree to which an analyte component is delayed as it is eluted through a given system (retentivity).

$$k' = (V_r - V_0)/V_0 = (t_r - t_0)/t_0$$

Where  $V_r$  = peak retention volume;  $V_0$  = column void volume  
 $t_r$  = peak retention time;  $t_0$  = peak void time

- The larger the  $k'$ , the later the analyte elutes after the void.

# Effect of $k'$ on Overall Resolution




- As  $k'$  grows larger, its effect reaches a limit at a value of about 10.
- Since  $k'$  depends on retention time, longer columns eventually have a diminished effect on resolution.

# Influencing the Capacity Factor $k'$


- Retentivity ( $k'$ ) decreases 2 - 3 fold for each 10% increase in mobile phase strength.
- Mobile Phase Strength -
  - As per the rule of thumb, altering the mobile phase strength also alters the retention of the analytes.
- Bonded Phase Functionality (Reverse Phase) -
  - As the bonded phase hydrophobicity increases (increasing alkyl chain length, etc.) so will the retention of the analytes.
- Temperature -
  - As temperature increases, the retention time decreases. This does not necessarily result in poorer separation because of the other factors in the resolution equation.

Which of these is easiest to change??


# Mobile Phase Strength vs. $k'$



# Temperature Effect on $k'$



# Temperature Effect on $k'$



# Summary of $k'$ Effects

- A larger value of  $k'$  means better resolution...to a certain extent ( $k' = 10$  maximum)
- Increasing the mobile phase strength decreases  $k'$
- Increasing the temperature decreases  $k'$ , but may not result in a “bad” separation based on the other factors affecting resolution



# Selectivity Factor, $\alpha$

- The selectivity or separation factor represents the ratio of any two adjacent  $k'$  values, thereby describing the relative separation of adjacent peaks. This relationship is expressed as:

$$\alpha = k'_b/k'_a$$

- If  $\alpha = 1$ , two components are perfectly overlapping
- For early eluting peaks you want  $\alpha$  to be large for good resolution.
- For later eluting peaks,  $\alpha$  can be smaller and still have acceptable separation.

# Effect of $\alpha$ on Overall Resolution

- Remember the resolution equation?

$$R_s = 1/4[(\alpha - 1) / \alpha] \sqrt{N} [k' / (1 + k')]$$

- Let's only look at the part involving  $\alpha$

$$R_s = 1/4[(\alpha - 1) / \alpha]$$


- And see how much resolution will improve with small changes in  $\alpha$

# Effect of $\alpha$ on Overall Resolution

$$R_s = 1 / 4 [(\alpha - 1) / \alpha]$$

- For an  $\alpha$  value of 1.1, the contribution of the selectivity term is
  - $(1.1 - 1) / 1.1 = 0.09$
- For an  $\alpha$  value of 1.4, the contribution of the selectivity term is
  - $(1.4 - 1) / 1.4 = 0.29$
- So...a very small change in  $\alpha$  leads to a more than THREE-FOLD increase in the contribution to resolution.

# Effect of $\alpha$ on Overall Resolution



- As  $\alpha$  grows larger, its effect reaches a limit at a value of about 5.
- Since  $\alpha$  depends on components' retention factor  $k'$ , longer columns eventually have a diminished effect on resolution.

# Influencing the Selectivity Factor $\alpha$

- Mobile Phase Type -

- The importance of the type of interactions between the mobile phase and analytes is critical to the optimization of the selectivity of a system.

- Column Type -

- The bonded phase functionality can be selected by its chemical nature to provide better selectivity in an analytical method.

- Temperature -

- Selective interactions between analyte molecules and the stationary phase may not become evident until a critical temperature is attained.

Which of these is easiest to change??

# Summary of $\alpha$ Effects


- Since  $\alpha$  is the ratio of two  $k'$  values, the same general statements apply:
  - Increasing the mobile phase strength decreases individual values of  $k'$ , but their ratio ( $\alpha$ ) may affect resolution
  - Increasing the temperature decreases individual values of  $k'$ , but their ratio ( $\alpha$ ) may significantly affect resolution.
- A small increase in  $\alpha$  leads to a large increase in resolution

# Column Efficiency, N

- The column efficiency is defined as the degree to which a column and/or other system components can physically and chemically affect the separation of analytes.
- As column efficiency increases, analyte components will elute in a smaller volume of the mobile phase, usually observed as narrower or “sharper” peak shapes.
- Column efficiency is generally expressed in terms of theoretical plate number.

# Calculation of Theoretical Plates

$$N = A(t_r / W)^2$$




| <u>W</u>      | <u>A</u> | <u>Method</u> | <u>Width measured at</u>                |
|---------------|----------|---------------|---|
| $W_i$         | 4        | Inflection    | Inflection point (60.7% of peak height) |
| $W_h$         | 5.54     | ½ Height      | 50% of peak height                      |
| $W_{3\sigma}$ | 9        | $3\sigma$     | 32.4% of peak height                    |
| $W_{4\sigma}$ | 16       | $4\sigma$     | 13.4% of peak height                    |
| $W_{5\sigma}$ | 25       | $5\sigma$     | 4.4% of peak height                     |
| $W_b$         | 16       | Tangent       | Baseline, following tangent drawing     |

Constants A are different at each peak width, assuming a perfect Gaussian shape.


Real-world peaks often have tailing, so widths measured at the lower part of the peak more accurately reflect the tailing when calculating N.



# Calculation of Efficiency, N



- Width measured at the baseline after tangent lines are drawn on the peak. Used when tailing is minimal.



- Width measured at 4.4% of peak height, no tangents drawn. Used when tailing is significant.

# Effect of N on Overall Resolution

- Do you STILL remember the resolution equation?

$$R_s = 1/4[(\alpha - 1)/\alpha]\sqrt{N}[k'/(1 + k')]$$

- Now let's look at the part involving N

$$R_s = 1/4\sqrt{N}$$

- And see how much resolution will improve with changes in N


# Effect of N on Overall Resolution

$$R_s = 1 / 4 \sqrt{N}$$

| <u>Plates</u> | <u><math>\sqrt{N}</math></u> | <u>Contribution</u> |
|---------------|------------------------------|---------------------|
| 5,000         | 70.7                         | - - - -             |
| 10,000        | 100                          | 41%                 |
| 20,000        | 141.4                        | 100%                |

- Since the contribution of N to resolution is a square root, doubling N from 5000 to 10,000 only increases the contribution to resolution by 41%.
- To double the effect on resolution coming from N, we have to increase the value of N by a factor of 4

# Effect of N on Overall Resolution



- Note that there is no flattening of the curve like with  $k'$  and  $\alpha$ .
- Resolution will continue to increase as theoretical plates increase.

# Influencing the Efficiency, N


- Particle Size and Size Distribution -
  - The smaller the particle size and the narrower the range of the particle size distribution, the more efficient the column.
- Packing Type -
  - Totally porous particles will also have greater efficiency than solid or pellicular-shaped packings, due to the additional surface area attributable to the pores.
- Mobile Phase Viscosity -
  - As mobile phase viscosity increases, molecular movement through the mobile phase is inhibited.
- Temperature -
  - For reverse phase chromatography, an increase in efficiency, N, may be realized as column temperature is increased.

# Effect of Particle Size on N

| Column Diameter (mm) | Column Length (cm) | Particle Size ( $\mu\text{m}$ ) | $4\sigma$ Peak Width ( $\mu\text{L}$ ) | Theoretical Plates per centimeter |
|----------------------|--------------------|---------------------------------|--|-----------------------------------|
| 10                   | 25                 | 10                              | 1118                                   | 333                               |
| 4.6                  | 25                 | 10                              | 237                                    | 333                               |
| 4.6                  | 25                 | 5                               | 167                                    | 667                               |
| 4.6                  | 10                 | 5                               | 106                                    | 667                               |
| 4.6                  | 10                 | 3                               | 82                                     | 1111                              |
| 4.6                  | 3                  | 3                               | 45                                     | 1111                              |
| 3                    | 10                 | 5                               | 45                                     | 667                               |
| 2                    | 25                 | 10                              | 45                                     | 333                               |
| 2                    | 25                 | 5                               | 32                                     | 667                               |
| 2                    | 10                 | 5                               | 20                                     | 667                               |
| 2                    | 10                 | 3                               | 15                                     | 1111                              |
| 1                    | 25                 | 10                              | 11                                     | 333                               |
| 1                    | 25                 | 5                               | 8                                      | 667                               |
| 1                    | 25                 | 3                               | 6                                      | 1111                              |
| 1                    | 10                 | 5                               | 5                                      | 667                               |
| 1                    | 10                 | 3                               | 4                                      | 1111                              |

Smaller particle sizes result in higher numbers of theoretical plates

# Summary: Review of Terms



# Summary: Relative Influence of All Factors on Resolution

| <u>Parameter Change</u> | <u>N</u>      | <u>k'</u> | <u><math>\alpha</math></u> | <u>R<sub>s</sub></u> |
|-------------------------|---------------|-----------|----------------------------|----------------------|
| <b>Standard</b>         | <b>10,000</b> | <b>2</b>  | <b>1.1</b>                 | <b>1.52</b>          |
| +10% N                  | 11,000        | 2         | 1.1                        | 1.59                 |
| -25% N                  | 7,500         | 2         | 1.1                        | 1.31                 |
| -50% N                  | 5,000         | 2         | 1.1                        | 1.07                 |
| -60% N                  | 4,000         | 2         | 1.1                        | 0.96                 |
| -75% N                  | 2,500         | 2         | 1.1                        | <u>0.76</u>          |
| +10% k'                 | 10,000        | 2.2       | 1.1                        | 1.56                 |
| +10% $\alpha$           | 10,000        | 2         | 1.2                        | <u>2.78</u>          |

Note that changing  $\alpha$  a very small amount has the biggest effect



# Summary: Review of Factors

| PARAMETER                | INFLUENCED BY                                     | TARGET VALUE   |
|--------------------------|---|--|
| Efficiency, $N$          | Column, System<br>Flowpath<br>Configuration       | Minimum of 400<br>Theoretical Plates per<br>centimeter |
| Capacity<br>Factor, $k'$ | Mobile Phase<br>Strength                          | 1.0 - 10   |
| Selectivity, $\alpha$    | Mobile Phase<br>Type,<br>Stationary Phase<br>Type | 1.1 - 2  |
| Resolution, $R_s$        | All of the Above                                  | 1.3 - 1.5 or Greater                                   |

# Questions About Columns?

## Next – HPLC System Components

# HPLC System Components

- Pumps
  - Micro to Analytical to Preparative Flow Rates
  - Isocratic and Gradient Configurations
- Degasser
  - How it Affects Pumping and Sample Injection
- Valves
  - Solvent Selection and Flow Selection

# HPLC System Components

- Sample Injection
  - Manual Injector or Autosampler
- Oven
  - How Temperature Affects Separation
  - Valves for Column Switching
- Detectors
  - UV-VIS
  - Diode Array
  - Fluorescence
  - Light Scattering
  - Refractive Index
  - Conductivity
  - Mass Spectrometer

# HPLC System Components

- Fraction Collector
  - Isolate Specific Sample Components
  - Purify Compounds for Multi-Step Synthesis
- Column
  - Types of Packing Material
  - Factors Affecting Separation
    - Particle Size and Column Length
    - Flow Rate and Temperature


# Hardware Components of an HPLC System

# HPLC Pumps – 2 Basic Types

- Tandem piston
  - Two pistons with different volumes (48 and 24  $\mu\text{L}$ )
  - During each stroke, 24  $\mu\text{L}$  of liquid is delivered
  - Best for higher analytical flow rates, up to 10 mL/min
  - Some pulsation is observed, and pulse dampeners are available
  - Not recommended for pulse-sensitive detectors like RID and CDD

# Tandem Piston Pump

Secondary Piston



Primary Piston




# HPLC Pumps – 2 Basic Types

## ■ Dual Piston

- Two pistons with equal volume (10  $\mu\text{L}$  each)
- During each stroke, 10  $\mu\text{L}$  is delivered
- Best for low flow rates ( $< 1 \text{ mL/min}$ )
- Little to NO pulsation, so it's ideal for pulse sensitive detectors like RID and CDD


# Dual Piston Pump



# Other Pump Components


- Check Valves

- Control liquid movement in and out of the pump head



# Other Pump Components

- Piston/plunger seal
  - Prevents solvent leakage out of pump head



Seal


- Inline filter
  - Removes solvent particulates

# HPLC Degassing

- Degassing removes dissolved air that interferes with check valve operation
- Helium sparge
  - Gas line from the tank directly in the solvent bottle
- Vacuum degassing
  - Sonicate before connecting to the system
  - Online with a degassing unit


# Valves Used With Pumps

- Solvent Selection – 2 Solvents Per Pump
  - Use for solvent switching




# Valves Used With Pumps


- Solvent Selection – 2 Solvents Per Pump
  - Use for pump loading of large sample volumes



Pump B – strong gradient solvent. Form the gradient with B.CONC command




Pump A – weak gradient solvent and sample loading



# Valves Used With Pumps

- Solvent Selection – 4 Solvents Per Pump
  - Use for low pressure gradient formation




Combine any proportion of A/B/C/D.

REQUIRES additional mixing before the injector.




# Valves Used With Pumps

- Solvent Selection – 4 Solvents Per Pump
  - Use for different gradients in method development



# Sample Injection – Manual


- Manual Injector with Syringe
  - Fixed loop of varying sizes (1 to 20 mL or more)
  - Fill with syringes of varying sizes
  - Can include a switch to start a data system



# Sample Injection – Automatic

- Fixed-Loop Autosampler
  - Loop is installed on the valve and can be changed for different injection volumes
  - External syringe draws sample and fills loop
- Advantages: low cost, rugged, few moving parts
- Disadvantages: Poor performance for low volume injections, higher carryover, always some sample loss

# Sample Injection – Fixed Loop




- External syringe draws sample, then fills the fixed-volume loop attached to the valve.


# Sample Injection – Automatic

- Needle-in-the-flowpath autosampler
  - Sample loop and needle are a single piece of tubing
  - Loop and needle are cleaned during the run
  - Metering pump draws sample very precisely
- Advantages: no sample loss, low carryover
- Disadvantages: higher cost, more delay volume for gradient

# Sample Injection to Flow Path





- Sample Loading



- Sample Injection – Everything drawn into the needle goes to the column.

# Rinsing After Injection



- Rinse liquid flows through ports 5 and 6 of the high pressure valve.
- Sample aspiration uses port 5.
- If air is present around port 5, injection reproducibility will be low.
- Rinse liquid **MUST** be degassed!


# HPLC Column Ovens

- Block heater with solvent preheater
  - Column is housed between 2 metal plates
  - Mobile phase is plumbed into the block for preheating
- Forced air
  - Column is in a large chamber with air circulation
  - Better temperature equilibration
  - Room for column switching valves



# Why Use a Column Oven?

- Retention times decrease, and higher flow rates are possible




# HPLC Detectors

- UV-VIS
- Diode Array
- Refractive Index
- Fluorescence
- Light Scattering
- Conductivity
- Mass Spectrometer

# HPLC Detectors – UV-VIS

- UV-VIS
  - Wavelength range 190-700 nm
  - D2 and W lamps
- Most common HPLC detector for a variety of samples
  - Proteins and peptides
  - Organic molecules
  - Pharmaceuticals
- Monitor 2 wavelengths at one time


# HPLC Detectors – UV-VIS



# HPLC Detectors – Diode Array

- Diode Array
  - Wavelength range 190-900 nm
  - D2 and W lamps
- Spectral information about sample
- Create compound libraries to identify unknowns
- Monitor an entire wavelength range at one time – up to 790 wavelengths vs. only 2 with a UV detector


# HPLC Detectors – Diode Array



# HPLC Detectors

- Refractive Index
  - For samples with little or no UV Absorption
    - Alcohols, sugars, saccharides, fatty acids, polymers
  - Best results when RI of samples is very different from RI of mobile phase
  - Flow cell is temperature controlled with a double insulated heating block.
- **REQUIRES** isocratic separations
- **REQUIRES** low pulsation pumps


# HPLC Detectors – RI Balance



- Fill sample and reference cell with mobile phase.




# HPLC Detectors – RI Analyze



- Mobile phase flows through sample side only.

# HPLC Detectors – RI Analyze




- As the refractive index changes, the image on the photodiode is deflected or “unbalanced”, and the difference in current to the photodiode is measured.

# HPLC Detectors

- Fluorescence
  - Xenon lamp for light source
  - Excitation wavelength range: 200-650 nm
  - Emission wavelength range: up to 900 nm depending on photomultiplier installed
- Used primarily for amino acid analysis
  - Derivatize samples before (pre-column) or after separation( post-column)


# HPLC Detectors - Fluorescence



# HPLC Detectors


- Evaporative Light Scattering (ELSD)
  - Also for low or no UV absorbing compounds
    - Sometimes called a “Universal” detector
  - Requires NO equilibration (unlike RID)
  - Can be used with gradients and volatile buffers (unlike RID)
  - Semi-volatile compounds can be detected at low temperatures

# ELSD Operation



# ELSD vs. Other Detectors

## Simultaneous ELSD-LT, RID and UV for Carbohydrates



- ELSD has higher sensitivity than UV and RID
- ELSD can be used with gradients, unlike RID

# HPLC Detectors

- Conductivity
  - Flow cell contains 2 electrodes
  - Measure ion amounts in sample
- **REQUIRES** low pulsation pumps
- Flow cell must be placed in a column oven



# HPLC Detectors - Conductivity

## ■ Conductivity

### ■ Use in Environmental and water testing

■  $\text{F}^-$ ,  $\text{Cl}^-$ ,  $\text{NO}_3^-$ ,  $\text{PO}_4^{3-}$ ,  $\text{SO}_4^{2-}$

■  $\text{Li}^+$ ,  $\text{Na}^+$ ,  $\text{K}^+$ ,  $\text{Mg}^{2+}$ ,  $\text{Cu}^{2+}$ , M-CN complexes

### ■ Determine organic acids in fruit juice

■ Oxalic, Maleic, Malic, Succinic, Citric

### ■ Analyze surfactants

■ Sulfonates, long/short chain ammonium


# HPLC Detectors

- Mass Spectrometer
  - Separate sample components as ions according to their mass to charge ( $m/z$ ) ratio
- Three stages to detection
  - Vaporization: liquid from HPLC column converted to an aerosol
  - Ionization: neutral molecules converted to charged species (either positive or negative)
  - Mass Analysis: filter ions by  $m/z$  ratio

# HPLC Detectors – Mass Spec

- Two Ionization Types
- APCI: Atmospheric Pressure Chemical Ionization
  - For molecules up to 1000 Da
  - Singly charged ions
  - Best for analysis of non-polar molecules
- ESI: Electrospray Ionization
  - Can be used for large biopolymers
  - Forms multiply charged ions
  - Best for the analysis of polar molecules, especially pharmaceutical products and proteins


# HPLC Detectors – Mass Spec



# HPLC System Components

## ■ Fraction Collector

- Purify raw materials or compounds from synthesis
- Collect by slope, level, time, volume
- Isolate single peaks per tube, or divide peaks into small “slices” for extra purity



# Questions About Hardware Components??

Next – HPLC System Types.





Now that we have hardware components and columns, what do we DO with them??

# HPLC System Types

- Isocratic system
  - Same mobile phase concentration throughout the separation
  - Use 1 pump and pre-mix solvents
  - Use 1 pump and a valve for 4 different solvents
  - Use 2 pumps and vary the amount coming from each pump

# Isocratic Separation

- 1 pump and premixing
- 4.6 mm ID Column, 1 mL/min, Changing MeOH % vs Water





# Isocratic Separation

- 1 pump with valve and premixing



A = 80% Methanol, 20% Water

B = 70% Methanol, 30% Water

C = 60% Methanol, 40% Water

D = 50% Methanol, 50% Water

# Isocratic Separation

- 1 pump with mixer – let the pump do the work!



Method 1: A.CONC = 20%, B.CONC = 80%


Method 2: A.CONC = 30%, B.CONC = 70%

Method 3: A.CONC = 40%, B.CONC = 60%

Method 4: A.CONC = 50%, B.CONC = 50%

# Low Pressure Gradient

- 1 Pump, solvents are mixed before the pump




- **REQUIRES** degassing

# HPLC System Types


- High Pressure Gradient
  - Multiple pumps are used with a mixer after the pumps
- Low Pressure Gradient
  - Solvents are mixed before the pump

# High Pressure Gradient

- Binary Gradient
- 2 Pumps and Mixer




- Ternary Gradient
- 3 Pumps and Mixer



# Low Pressure Gradient

- 1 Pump, solvents are mixed before the pump



- **REQUIRES** degassing

**Questions About System Types?**

**Next: Troubleshooting and How  
to Take Care of Your Column and  
HPLC System**

# HPLC Troubleshooting

- Pressure: too much or too little
- Leaks: pump, autosampler, detector
- Reproducibility: pump, autosampler
  
- Column Care: Flushing and equilibration



# Pump Troubleshooting

- No pressure, or fluctuating pressure
  - Pump may not be completely full of liquid – check solvent inlet line
  - Air in check valve – always degas mobile phase!
  - “Stuck” check valve – the pump may have been idle for too long and solvent has dried inside the check valve.

Poor quality solvent: may contain resins that coat the ball inside the check valve, and that film won't let the ball seat properly

# Pump Troubleshooting

- High Pressure
  - Outlet frit may be blocked with particles from mobile phase or seal material
- Leaks
  - Damage to seal and/or plunger due to several factors
    - Misaligned plunger
    - Solvent incompatibility with seal material
    - Salt crystal buildup from buffers – use a rinse kit!

# Pump Troubleshooting

- Retention Time Reproducibility
  - For a dual piston pump, only one side may be filled with liquid – check solvent inlet lines
  - Temperature change (may not be the pump's fault)
    - A 1° shift in temperature can result in a 1-2% shift in retention time
    - Avoid drafty locations in the lab
    - Use a column oven when possible

# Autosampler Troubleshooting

## ■ High Pressure

- Particulates from mobile phase, sample, pump may be trapped in the inlet tubing or valve
  - Filter mobile phase AND sample when possible

## ■ Leaks

- Fittings may be loose on the valve
  - Tighten fittings properly and don't exceed the pressure limit of the autosampler

# Autosampler Troubleshooting

- Area % Reproducibility
  - Always degas rinse phase, and use some volume of liquid for rinsing to keep all flow paths in the valves full of liquid
  - Make sure the needle stroke is deep enough to draw sample from the vial
  - Check for leaks on the valve fittings, and the connection to the column inlet

# Detector Troubleshooting

## ■ Spiky Baseline

- Air bubble in flow cell – degas mobile phase!
  - Put some restriction on the cell outlet, but not too much!  
Tubing with 0.005” i.d. is fine.

## ■ Leaks

- Cracked flow cell
  - Don't exceed the pressure limit of the cell
- Poor tubing connections
  - Use the proper fittings and tighten appropriately

# Column Care

- Follow MFR's recommendations for solvent compatibility, flow rate, and pressure limits
- Filter samples when possible
  - Particulates will build up on the inlet frit over time
- Use care when reversing column flow
  - Connect the outlet to waste, NOT inline with the detector to prevent further contamination
- Store columns in recommended solvents

# Troubleshooting Summary

- Throw away bad parts and columns.
- Leaks do not fix themselves.
- If it doesn't pass, you must degas.



Questions About  
Troubleshooting?

Tomorrow: Application-Specific  
Systems, Software, and  
Prominence Demonstration

# HPLC Applied Systems

## ■ Protein Separations

- Column selection is important: reversed phase C-18, ion exchange most common
- Buffered mobile phases often used so a rinse kit for the pumps is recommended
- Inert (PEEK) pump and autosampler may be necessary
- UV or Diode Array detection
- Fraction collection for isolation and purification


# HPLC Applied Systems

## ■ Proteomics


- Very small sample amounts with many components
- Use 2-dimensional chromatography
  - Elute portions of sample onto a trap column with a salt gradient
  - Desalt the trap then transfer sample to reversed phase column
  - Elute with a reversed phase gradient

# 2-Dimensional HPLC

Load sample to SCX Column and elute portion to Trap




# 2-Dimensional HPLC



# 2-Dimensional HPLC

Flush Trap to RP Column and elute to detector



# HPLC Applied Systems

- Amino Acid Analysis
  - Column selection is important: C-18 is very common
  - Any pumps, autosampler, oven
  - Pre- or post column derivatization (OPA)
    - Autosampler can do pre-column reactions
    - Additional pump for post-column reagent addition
  - Fluorescence detection most common


# HPLC Applied Systems

- Food and Beverage Industry
  - Many isocratic methods
  - C18 columns, ion exchange columns
  - Any pumps, autosampler, oven
  - Traditional methods use UV, RID
  - Perfect opportunity for ELSD: App. notes on
    - Chili peppers
    - Wine
    - Sugar alcohols
    - Cereal



# ELSD for Food and Beverage

## Mono-, Di- and Oligosaccharide Standards




1. Glycerol
2. Arabinose
3. Fructose
4. Glucose
5. Sucrose
6. Maltose
7. Maltotriose
8. Maltotetraose
9. Maltopentaose
10. Maltohexaose
11. Maltoheptaose

**Column:** Asahipak NH<sub>2</sub>-P50, 5 $\mu$ m, 250x4.6mm  
**Mobile Phase:** A:Acetonitrile B: 0.0004N NH<sub>4</sub>OH  
**Gradient:** (Time, %B)(0,15)(60,65)  
**Flowrate:** 1.0mL/min  
**Col. Temp:** 30°C  
**Detector:** Shimadzu ELSD-LT (Gain 5; T 40°C; P 250kPa)

# ELSD for Food and Beverage


## Fructo-oligosaccharides in Reconstituted, Minced Onions



Column: Asahipak NH<sub>2</sub>-P50, 5 $\mu$ m, 250x4.6mm  
Mobile Phase: A:Acetonitrile B: 0.0004N NH<sub>4</sub>OH  
Gradient: (Time, %B)(0,15)(60,65)  
Flowrate: 1.0mL/min  
Col. Temp: 30°C  
Detector: Shimadzu ELSD-LT (Gain 6; T 40°C; P 250kPa)

# ELSD for Food and Beverage

## Malto-oligosaccharide Profile of Domestic Ale



Column: Asahipak NH<sub>2</sub>-P50, 5 $\mu$ m, 250x4.6mm  
Mobile Phase: A:Acetonitrile B: 0.0004N NH<sub>4</sub>OH  
Gradient: (Time, %B)(0,15)(60,65)  
Flowrate: 1.0mL/min  
Col. Temp: 30°C  
Detector: Shimadzu ELSD-LT (Gain 8; T 40°C; P 250kPa)

# HPLC Applied Systems

- Nutraceutical: \$46.7 BILLION In 2002, predicted to grow almost 10% each year\*.
- Watch for these keywords
  - Functional foods/beverages
  - Fortified
  - Energy/nutrition
  - Health-promoting
  - Natural/Herbal
  - Vitamin/Mineral/Supplement

\* <http://www.bccresearch.com/editors/RGA-085R.html>


# HPLC Applied Systems

- Nutraceutical system configurations
  - Similar to Food and Beverage
  - Promote ELSD since many compounds have low (or no!) UV absorbance
- There are many application notes available for nutraceutical samples
  - White Willow Bark
  - Black Cohosh
  - Milk Thistle


# ELSD for Nutraceutical

## Siberian Ginseng Extract

Brand A Label Claim:  
No Excipients Present



Brand B Label Claim:  
Numerous Excipients




0 4 8 12 16 20 24 28 32 36 40 44 48 52  
Minutes

Column: Asahipak NH<sub>2</sub>-P50, 5 $\mu$ m, 250x4.6mm  
Mobile Phase: A:Acetonitrile B: 0.0004N NH<sub>4</sub>OH  
Gradient: (Time, %B)(0,15)(60,65)  
Flowrate: 1.0mL/min  
Col. Temp: 30°C  
Detector: Shimadzu ELSD-LT (Gain 5; T 40°C; P 250kPa)


# ELSD for Nutraceutical

Saw Palmetto Oil



1. LGG
2. LGL
3. LLL
4. OGL
5. PGL
6. OLL
7. OGO

Evening Primrose Oil



**Column:** Shimadzu Premier C18, 5 $\mu$ m, 150x4.6mm  
**Mobile Phase:** A: Acetonitrile B: Dichloromethane  
**Gradient:** (Time, %B)(0,15)(20,30)(40,70)  
**Flowrate:** 1.0mL/min  
**Col. Temp:** 30°C  
**Detector:** Shimadzu ELSD-LT

# HPLC Applied Systems

- Ion Chromatography
  - Column selection is most important
  - Low pulsation pumps and any autosampler
  - UV or Conductivity detector
    - Ion chromatography applications data book
    - Suppressed or non-suppressed detection
      - Metrohm-Peak Model 833
      - Alltech Model 640 or 641



# Ion Chromatography Applications


- Inorganic Anions – tap water
  - $F^-$ ,  $Cl^-$ ,  $NO_3^-$ ,  $PO_4^{3-}$ ,  $SO_4^{2-}$
- Cations and Transition Metals – tap water
  - $Li^+$ ,  $Na^+$ ,  $K^+$ ,  $Mg^{2+}$ ,  $Cu^{2+}$ , M-CN complexes
- Organic Acids – fruit juice
  - Oxalic, Maleic, Malic, Succinic, Citric
- Surfactants – soaps and detergents
  - Sulfonates, long/short chain ammonium

# Ion Chromatography Columns

- Alltech
- Phenomenex
- Dionex
  
- Silica and polystyrene-based with specific functional groups

# Ion Chromatography Applications

## Common Cations




|              | (ppm) |
|--------------|-------|
| 1. Potassium | 2.5   |
| 2. Magnesium | 2     |
| 3. Calcium   | 2     |
| 4. Ammonium  | 1.5   |
| 5. Sodium    | 1.5   |
| 6. Lithium   | 0.2   |

**Column:** ShimPak IC-C3, 5 $\mu$ m, 150x4.6mm  
**Mobile Phase:** 2.5mM oxalic acid  
**Flowrate:** 1.5mL/min  
**Col. Temp.:** 40 °C  
**Cell Temp.:** 43 °C  
**Inj. Vol.:** 30 $\mu$ L  
**Detector:** Shimadzu CDD-10AVP non-suppressed  
(Gain 2; Polarity -1; Response 4)

# Ion Chromatography Applications

## Common Anions




|             | A(ppm) | B(ppm) |
|-------------|--------|--------|
| 1. Fluoride | 25     | 0.6    |
| 2. Chloride | 50     | 1.3    |
| 3. Nitrite  | 50     | 1.3    |
| 4. Bromide  | 50     | 1.3    |
| 5. Nitrate  | 50     | 1.3    |
| 6. Sulfate  | 50     | 1.3    |

**Column:** ShimPak IC-A3, 5 $\mu$ m, 150x4.6mm  
**Mobile Phase:** 2mM phthalic acid @pH 4.2 with LiOH  
**Flowrate:** 1.5mL/min  
**Col. Temp.:** 37 °C  
**Cell Temp.:** 40 °C  
**Inj. Vol.:** 10 $\mu$ L  
**Detector:** Shimadzu CDD-10AVP non-suppressed  
(Gain 2; Polarity 1; Response 4)


# Research and Method Development

- Typically, more “advanced” systems use multiple detectors and valves for column and solvent switching



# Research and Method Development

- Some “advanced” systems will include a high capacity autosampler and a mass spectrometer



# Application Questions?

## Next: Software Demonstration and Prominence Hardware

# Prominence Overview

- System Controller
- Pump and Degasser
- Autosampler and Rack Changer
- Column Oven and Valves
- UV and Diode Array Detectors



# CBM System Controller

- Web-based control
  - Connect to lab network or directly to computer
- Methods stored in CBM or connected computer
- Controls all components that have a fiber optic cable
  - 10A and VP Series



# Standard Pump

- LC-20AT
  - 1  $\mu\text{L}$  to 10 mL/minute
  - LPGE valve can be installed in the pump
  - Reduced delay volume
  - Sapphire piston and GFP seal
  - Floating piston design




# Micro-Flow Pump

- LC-20AD
  - 0.1  $\mu\text{L}$  to 10.0 mL/min
  - 10  $\mu\text{L}$  pistons for no pulsation
    - RID, ECD, CDD
  - Sapphire piston and GFP seal
  - Ideal for low flow rate and LCMS applications



# Binary Pump

- LC-20AB
  - 2 LC-20AD in 1 box
  - Binary, space saving configuration
  - 0.1 to 10.0 mL/min
  - For gradient flow rate  $> 0.4$  mL/minute



# DGU-20A3 and A5 Degasser

- Vacuum degasser
  - Internal volume of  $< 400 \mu\text{L}$
  - Teflon AF membrane for efficient  $\text{O}_2$  removal
  - Plug into pump for power and control
  - External power supply available

# Autosampler

- Two Models:
  - SIL-20A
  - SIL-20AC: 4-40C temp. control
- Enhanced Carryover Performance
- Faster Cycle Time
- Optional Active Rinsing
- Optional Rack Changer






# Rack Changer

- Two Models
  - A; ambient or C; 4-40° C,  $\pm 6^\circ$  temp. control
- 12 x 96 well MTP racks (reg. or deep well) in 4 stacks
  - Mix and match plate type between stacks
  - ~90 seconds to change plates.



# Column Oven

- Forced air heating and cooling
  - CTO-20A: ambient – 85
  - CTO-20AC: (ambient -15) – 85
- Higher T.MAX for polymer and carbohydrate applications
- Linear temperature programming possible
- Integrated valve controller
- Space inside for 2 switching valves





# Switching Valves

## ■ FCV-20AH<sub>2</sub>

- 2 Position 6 port High Pressure valve
- Column Switching
- Standalone control possible (front panel or Event) OR install in CTO-20A/AC

## ■ FCV-20AH<sub>6</sub>


- 6 Position 7 port High Pressure valve
- Column Selection
- Standalone control possible (front panel of Event) OR install in CTO-20A/AC

# UV Detector

- Extended wavelength range (190-700 nm)
- Improved Noise and Drift Specs
- Temp Controlled Flow Cell
- 2.5 AU Linear Range
- Included Hg lamp for wavelength accuracy



# Thermostatted Flow Cell



Effect of Temperature Controlled Cell — Room Temperature raised from 20-30C

# Diode Array Detector

- World's lowest noise PDA
- World's best linearity -  $\geq 2.0$  AU
- Temperature Controlled Flow Cell
- Variable Slit Width
  - 8 nm (better S/N) and 1.2 nm (better resolution)
- 4 Channel Analog Board is STD
- Ethernet Communication

# LC-2010 Integrated HPLC System

- Fully integrated HPLC system ideal for:
  - QA/QC environment
  - High-throughput applications
  - University teaching laboratories
  
- Standalone or software controlled
  - Easy to navigate control screens
  - GUI with “Wizard” assistance
  - Standard or “simple” mode

# LC-2010HT Features



- Dynamic inlet valve
  - Quaternary gradient unit
- High speed autosampler
  - 4-40 C temperature control
- Column heater
- 2.5 AU detector linearity
  - Thermostatted flow cell
- Automatic power, system prep, and validation functions

# LC-2010HT Pumping System

- 5-channel degassing unit
  - 4 mL/line for solvents A-D, 2 mL/line for SIL
- Dynamic Inlet Valve
  - Electronic check valve to keep prime and minimize air bubbles
- 4 solvent proportioning valve (FCV-10ALvp style)
  - Gradient accuracy of  $\pm 0.5\%$
- Manual or automatic priming



# LC-2010HT Pump Performance

- Units are pre-plumbed; users only add a column
- Instrument-to-instrument uniformity
  - 7 instruments, same column and paraben test mixture

| Mean retention time, 6 reps |      | Methyl   | Ethyl    | Propyl  | Butyl    |
|-----------------------------|------|----------|----------|---------|----------|
| S/N 005                     |      | 1.693    | 2.217    | 3.245   | 4.457    |
| S/N 051                     |      | 1.680    | 2.192    | 3.197   | 4.402    |
| S/N 054                     |      | 1.677    | 2.185    | 3.177   | 4.387    |
| S/N 056                     |      | 1.698    | 2.222    | 3.242   | 4.450    |
| S/N 058                     |      | 1.690    | 2.208    | 3.222   | 4.415    |
| S/N 060                     |      | 1.687    | 2.205    | 3.220   | 4.422    |
| S/N 062                     |      | 1.670    | 2.173    | 3.155   | 4.363    |
|                             |      |          |          |         |          |
|                             | %RSD | 0.581481 | 0.804374 | 1.04587 | 0.756171 |



# LC-2010HT Autosampler

- High Capacity
  - 350 1 mL vials, 210 2 mL vials (LC-2010A), 4 microtiter plates (96 and 384 well; Std or Deep-well)
- Fast injection
  - 15 second injection, ~30 second cycle time
- Reproducibility < 0.3% RSD specification
  - Typical value: ~0.10%
- Low carryover: < 0.01% (naphthalene analysis)
  - **NEW** Pt coated needle, PEEK rotor and PEEK needle seal to further reduce carryover

# LC-2010HT Autosampler Performance

## ■ Injection Reproducibility

- Method: Isocratic premixed 60:40 MeOH:H<sub>2</sub>O
- Sample: Paraben test mix; 1, 5, 10, 25, and 50 µL injections, 10 reps each

| LC-2010A       | 1 µL %RSD | 5 µL %RSD | 10 µL %RSD | 25 µL %RSD | 50 µL %RSD |
|----------------|-----------|-----------|------------|------------|------------|
| methyl paraben | 0.295     | 0.0549    | 0.0393     | 0.0685     | 0.0425     |
| ethyl paraben  | 0.228     | 0.0705    | 0.0385     | 0.0370     | 0.0560     |
| propyl paraben | 0.327     | 0.0533    | 0.0509     | 0.0233     | 0.0463     |
| butyl paraben  | 0.285     | 0.0773    | 0.0336     | 0.0376     | 0.0439     |

| LC-2010C       | 1 µL %RSD | 5 µL %RSD | 10 µL %RSD | 25 µL %RSD | 50 µL %RSD |
|----------------|-----------|-----------|------------|------------|------------|
| methyl paraben | 0.283     | 0.0562    | 0.0392     | 0.0223     | 0.0515     |
| ethyl paraben  | 0.265     | 0.0533    | 0.0335     | 0.0325     | 0.0473     |
| propyl paraben | 0.246     | 0.0511    | 0.0511     | 0.0427     | 0.0405     |
| butyl paraben  | 0.265     | 0.0310    | 0.0204     | 0.0206     | 0.1210     |

# LC-2010HT Autosampler Performance

- Injector cycle time is crucial for high-throughput and mass spec. applications
- The LC-2010HT can inject in ~15 seconds
- Actual time, from pressing RUN to injection

|          | Rep 1 | Rep 2 | Rep 3 | Rep 4 | Rep 5 | Rep 6 |
|----------|-------|-------|-------|-------|-------|-------|
| LC-2010A | 15.58 | 15.49 | 15.24 | 15.52 | 15.43 | 15.64 |

# LC-2010HT Autosampler Performance

- Injection linearity
  - Paraben test mix: 1, 5, 10, 25, 50  $\mu\text{L}$  injections
  - 10 repetitions per level

| LC-2010A       | $R^2$    |
|----------------|----------|
| methyl paraben | 0.999976 |
| ethyl paraben  | 0.999982 |
| propyl paraben | 0.999987 |
| butyl paraben  | 0.999989 |


| LC-2010C       | $R^2$    |
|----------------|----------|
| methyl paraben | 0.999997 |
| ethyl paraben  | 0.999994 |
| propyl paraben | 0.999995 |
| butyl paraben  | 0.999994 |

# LC-2010HT Column Oven

- Block style that heats and cools column
  - Setting range of (Ambient - 15) to 60 C
- Adjustable aluminum blocks for extra contact points with column
- Solvent preheater: 4 or 9  $\mu\text{L}$
- Mixer in direct contact with heating block
  - Mixer volume is 240  $\mu\text{L}$

# LC-2010HT Detector

- 2.5 AU linearity spec
- Built in Hg lamp for wavelength calibration
- Thermostatted flow cell: 40 and 50 C settings
  - Prevents change in absorbance due to refractive index change with temperature variations



# LC-2010HT Detector Performance

- Linear to 2.5 AU
- Prednisone: 5 concentration levels
  - 10  $\mu$ L injections, 5 reps at each level
  - 60:40 MeOH:H<sub>2</sub>O, 4.6x100mm C18 column

| <b>LC-2010C S/N 002</b> |             |          |                  |
|-------------------------|-------------|----------|------------------|
| Level 1                 | 0.045 mg/mL | 144 mAU  |                  |
| Level 2                 | 0.090 mg/mL | 294 mAU  |                  |
| Level 3                 | 0.180 mg/mL | 600 mAU  |                  |
| Level 4                 | 0.360 mg/mL | 1215 mAU |                  |
| Level 5                 | 0.720 mg/mL | 2550 mAU | $R^2 = 0.999800$ |
| <b>LC-2010C S/N 003</b> |             |          |                  |
| Level 1                 | 0.045 mg/mL | 141 mAU  |                  |
| Level 2                 | 0.090 mg/mL | 288 mAU  |                  |
| Level 3                 | 0.180 mg/mL | 585 mAU  |                  |
| Level 4                 | 0.360 mg/mL | 1190 mAU |                  |
| Level 5                 | 0.720 mg/mL | 2530 mAU | $R^2 = 0.999688$ |

# Additional LC-2010HT Features

- Automatic power on/off
- System Prep – for running samples with different solvents
- Automatic system validation
- Individual component validation
- Status screen – mobile phase calculation
- Maintenance/parts replacement