

## Suggestions for JAG Group laboratory Notebooks

The following is a list of things that all laboratory notebooks should contain. It is by no means a comprehensive list and various experiments may require different forms of record keeping. This list is mainly meant for new initiates to learn what kinds of things are expected from them in their lab notebooks.

General laboratory notebook guidelines include:

1. A Table of Contents
  - a. This should include the experiment titles and the pages where they can be found.
2. A unified labeling scheme
  - a. A unified labeling scheme for experiments, data and compounds is useful for figuring out where things came from in relation to the notebook. The commonly accepted format is INITIALS-BOOK #-EXPERIMENT #.
  - b. For example, data labeled ALE-2-034, or ALE2034 would represent data for Alex Estrada, book 2, page 034.
3. All experiments should have page numbers along with descriptive/useful titles.
4. Experiments should include dates for ALL days an entry was made
5. Drawn reactions schemes should be as detailed as possible
6. References for experiments should be included.
7. There should be some reasoning or background information for the experiment.
8. Reagents should be in a table describing the amounts needed/used, their source, as well as any relevant safety issues.
  - a. If reagents have been prepared, there should be either an explanation or a reference to an explanation of how they were prepared.
9. Experimental sections should be as detailed as possible. They should contain all observations, measurements and procedures.

- a. This includes things like masses used, orders of addition, solution concentrations. Everything is important!
  - b. Don't forget about more subtle things like column dimensions and crystallization conditions.
  - c. Drawings of special setups or glassware are especially useful for reproducing an experiment.
  - d. TLC plates may also be put into this section. It may be useful to cover them completely in packing tape to prevent fading.
  - e. Experimental sections are also an appropriate place to mention safety hazards including, but not limited to, dangerous synthetic procedures or flammable/toxic compounds.
10. Data sections where all data is presented and worked up (if applicable) should also be present.
- a. For most synthetic experiments, data should include pertinent NMR signals, IR bands, etc.
  - b. Instruments and their locations should also be noted if they do not belong to the group.
  - c. Instrument parameters should also be noted here if different from typical values.
11. If data cannot be put in the notebook, there must be some clear indication as to where it can be found.
12. Some commentary or conclusions indicating what the results meant.
13. All notebooks should be written clearly and neatly. This includes having good penmanship as well as correct grammar.

BOOK # - PAGE #

EXPERIMENT TITLE AND CODE

REACTION  
SCHEME

Reference(s)

PURPOSE

PRELIMINARY CALCULATIONS

REAGENTS	SOURCE	AMOUNT	OTHER DETAILS

REAGENTS	SAFETY ISSUES

EXPERIMENTAL, DATA AND  
OBSERVATIONS

SIGNATURE



DATE

CONT. FROM PG.

BOOK # - PAGE #


EXPERIMENT TITLE AND CODE

EXPERIMENTAL, DATA AND OBSERVATIONS



DATA WORK-UP  
CHARTS, GRAPHS, TABLES

	A	B
1	###	###
2	###	###
3	###	###
4	###	###




SIGNATURE

DATE

TITLE Borane Protection of Dimethylamine

PROJECT ADL-2-121

Continued from page X



Ref: M. Van der Schelde et al. / Tetrahedron Lett (2009) 50:640-645

Purpose: Needed for the synthesis of "gearbox ligand"


SAFETY

- HPPH<sub>2</sub> - Causes skin & eye irritation (wear goggles & gloves)  
 - is pyrophoric & self-heating (avoid contact w/ air)
- H<sub>2</sub>B-SMe<sub>2</sub> - Flammable (avoid open flames), HAS acute oral & dermal toxicity so avoid contact  
 - releases flammable gases when contact H<sub>2</sub>O
- HPPH<sub>2</sub>-BH<sub>3</sub> - releases flammable gases in contact w/ H<sub>2</sub>O  
 - Store in a dry place but is air stable

Reagent	Source	MW	M	V	mass (g)	mmol
HPPH <sub>2</sub>	Strem	186.19	-	-	0.4137	2.222
H <sub>2</sub> B-SMe <sub>2</sub> (anHF)	Acros Org.	N/A	2	11.1mc	-	22.2

PROCEDURE

- Put hot 100ml schlenk tube, equipped w/ stirbar in g. box attached to org under vac. Then move into g. box
- ing. bar weigh 0.4137g of dimethylamine. Dissolve line of PCM + transfer to ~~flash~~ to flash. Remove flash from g. box.
- Put flash on line + flush with 3x before opening to N<sub>2</sub>
- Put bottle of H<sub>2</sub>B-SMe<sub>2</sub> online. Connect to N<sub>2</sub> + also insert a purge needle
- Flush a 24mc syringe 3x w/ N<sub>2</sub>.
- via syringe take 11.1mc of BH<sub>3</sub>SMe<sub>2</sub> soln + transfer purge needle to 100ml schlenk tube.
- Let stir RT. Rxn start @ 10:45 AM
- solution is clear + colorless



Continued to page 21

SIGNATURE

DATE  
 2/9/15

DISCLOSED AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued from page 121

2-10-15

- degassed anhydrous saturated NH<sub>4</sub>Cl soln for 1 hr.  
 - @ 10:30 led rxn mixture & slowly added 1 cm<sup>3</sup> of Nthyl Cl solution via syringe. Placed a needle in the septa in case too much pressure builds up (gas is produced)

- frequently produces white solid + lots of bubbly

→ this process is quiet exothermic. Take caution

→ finished @ 12:40 pm. removed vial both. 2151

- @ 2:00 pm dumped rxn mixture in 100 mL D<sub>2</sub>O

→ extracted w/ 3x 90 mL DCM | 750 mL sep funnel.

→ Then washed the organic layer (DCM) w/ 100 mL total saturated NaHCO<sub>3</sub> aq solution.

(Did this by dividing the collected layer into 2 portions.

Rinsed each portion w/ some of the NaHCO<sub>3</sub>.

→ Dried the DM solution w/ MgSO<sub>4</sub>

→ Filtered. used ~25 mL to rinse the MgSO<sub>4</sub>

→ reduced volume via volume removal solvent. (Carbon)


→ Filtered soln through 1 in pad of silica. rinse of the pad w/ 200 mL DCM. Transferred most of the soln to a 250 mL RB lit all wouldn't fit but saved w/ water fit.)

→ degassed soln w/ N<sub>2</sub>. ~~the soln was low enough added~~

~~the solution that didn't fit back. Let continue to evaporate~~

→ let soln that wouldn't fit evaporate in open air. Label ADL-2-121 B

bubbled N<sub>2</sub> through the soln (This will be a test to see how air sensitive it is)



2-22-15

- Solvent from RB has finally evaporated.

- redissolved in minimal DCM + transfer to weighed vial (not labeled) mass = 13.2415g

- removed DCM via rotary evaporator

(H<sub>2</sub>O bath @ ~30°C - 35°C)

Continued to page

SIGNATURE

DATE


DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE Synthesis of trans-Mo(CO)<sub>4</sub>(PPh<sub>2</sub>Cl)<sub>2</sub> PROJECT ABC-2-123

Continued from page 1



Reference: J. R. Martin et al / Journal of Organometallic Chemistry 751 (2014) 695-702

Purpose: To synthesize trans-Mo(CO)<sub>4</sub>(PPh<sub>2</sub>Cl)<sub>2</sub> as a Mo(CO)<sub>4</sub> source.

SAFETY

- Mo(CO)<sub>6</sub> - Poisonous (swallowed), in contact (skin or inhaled)
  - remove gloves, lab coat, do not touch
  - considered to have acute toxicity (cat. 2)
- PPh<sub>2</sub>Cl - Irritant, can cause severe skin + eye damage (handled in box so should not be a problem)

Reagents	Source	mw	mass (g)	mmol
Mo(CO) <sub>6</sub> (white powder)	STEM	264.00	<del>3.3221</del>	12.58
PPh <sub>2</sub> Cl (colorless liquid)	TGI	220.03	3.9434	17.87

Procedure

- Placed a hot 250 mL 2N <sup>N<sub>2</sub></sup> adapter on line. Flask was equipped w/ hot stir bar, reflux condenser + N<sub>2</sub> adapter @ top of reflux condenser. Let cool under vac thru rot valve system under N<sub>2</sub>.

- Weigh 3.3221g of Mo(CO)<sub>6</sub> + place in the flask. Flush flask 3x w/ N<sub>2</sub> + then draw under vac.

- Tear down septa + lines for flask under vac into the <sup>glass</sup> ~~glass~~ flask.

- In glovebox weigh 3.9434g of PPh<sub>2</sub>Cl in a vial. Transfer to 250 mL RB. Rinse vial w/ 4 mL Tol + transfer this to RB as well. Rinse flask w/ 4 mL Tol + transfer this to RB as well.

- Rinse flask w/ 4 mL Tol via syringe tech. from storage flask + into RB.

RB starts 3:00 PM

heat @ 110°C - solid dissolves almost immediately solution

LIQ + colorless  
Continued to page 124

SIGNATURE  
HOM Cebral

DATE  
2-9-15

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued from page 123

- Increase temp to 129°C - reflux does not occur

~~temp~~

- Increase temp to 125°C

- After heating for 10m rxn is now slightly yellow. reflux still has not begun.

- Increase temp to 130°C

- 25min in rxn is bright yellow still not refluxing though.

- increase temp to 135°C. still not refluxing. increase to 140°C

- Oranage by 35min in.

- hose to N<sub>2</sub> piped off. Wiggling it it back on & return temp to 135°C.

- ~~to rxn~~ boiling profusely but does not appear to be "refluxing"

- Slow reflux is 3:45 (1 drop/2min)

- rxn is red

- @ 4: S:O<sub>3</sub> rxn is dark red black + orange. realized tel was condensing in N<sub>2</sub> hose.

- attached another hose to 1000' condenser, flush 3x w/ N<sub>2</sub>, open to N<sub>2</sub>

- Close side adapt. & remove hose.

- let continue to reflux.

2-10-15

- in the morning mixture is black and is still refluxing slowly (still 135°C.)

- the oil bath has turned yellow. It does not appear to be a Mg fault though (no trace of metal left hose or condenser but, not sure what the contaminant is)


- at 8:30pm remove oil bath, let can cool, then remove solvent via vacuum pump. ~~For rxn under N<sub>2</sub> bar~~

- This produces a thick black residue. R445 stem under N<sub>2</sub> before leaving.

2-11-15

- Attempted to purify mixture w/ column chrom. used column w/ 10% Sy of silica + 3:1 DCM:hex solvent

- However nothing really separated. It pulled off as a big black smudge



SIGNATURE

*Am. Control*

DATE

2-11-15

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued to page 127